

Cassy

Closed Suction System

Intended use:

Closed Suction System Cassy is intended for removal of accumulated secretion in the lower respiratory tract in patients receiving ventilator treatment.

The product is intended for single use only.



Product range:

	Catheter CH	Outer diameter, mm	Length of suction catheter** cm			
			25	35	45	55
			REF*	REF*	REF*	REF*
Closed Suction System Cassy	06	2,0		3600118x	3600218x	
	08	2,7		3600418x	3600518x	3600618x
	10	3,3	3604518x	3600718x	3600818x	3600918x
	12	4,0	3604618x	3601018x	3601118x	3601218x
	14	4,7	3604718x	3601318x	3601418x	3601518x
	16	5,3	3604818x	3601618x	3601718x	3601818x
Closed Suction System Cassy with corrugated tube	10	3,3	3605118x	3601918x	3602018x	3602118x
	12	4,0	3605218x	3602218x	3602318x	3602418x
	14	4,7	3605318x	3602518x	3602618x	3602718x
	16	5,3	3605418x	3602818x	3602918x	3603018x

* Here noted as x, reference codes contain additional digits (1-5) representing a language code.

** Effective (usable) length of the suction catheter

Storage and handling:

